

Product Information

KinaseSTAR JNK Activity Screening Kit

I. Kit Contents:

Components	K2078-40 40 assays	Cap Color	Part Number
Kinase Extraction Buffer	80 ml	NM	K2078-C-1
c-Jun (1-79) Fusion Protein Beads	800 μ l	Red	K2078-C-2
Kinase Assay Buffer	25 ml	WM	K2078-C-3
ATP (10 mM)	50 μ l	Yellow	K2078-C-4
Phospho-cJun Specific Antibody	50 μ l	Green	K2078-C-5

II. Introduction:

c-Jun N-terminal kinase (JNK) is a member of the mitogen-activated protein (MAP) kinase family identified in mammals. Activation of JNK plays an critical role in neuronal apoptosis and other physiological and pathological processes.

The KinaseSTAR JNK Activity Screening Kit provides a simple and convenient way for screening of JNK activity based on Western Blot method. The assay uses an N-terminal c-Jun (1 - 79) fusion protein bound to glutathion sepharose beads to selectively “pull down” JNK from cell or tissue lysates. After washing to remove nonspecifically bound proteins, the kinase reaction is conducted in the presence of cold ATP. The phosphorylation of c-Jun is detected by Western blot method using a phospho-c-Jun specific antibody. The kit is suited for both research and drug discovery.

III. JNK Activity Immunoassay Protocol:

A. Preparation of Cell Lysate:

1. Activate cells by desired methods. Concurrently incubate a negative control culture without activation. To generate a positive control, cells can be treated with 1 μ g/ml of Anisomycin for 1 hr, before harvested.
2. Pellet cells (2-10 millions/assay) and wash once in 1X ice-cold PBS.
3. Lyse cells in 200 μ l ice-cold JNK Extraction Buffer. Incubate on ice for 5 min.
4. Pellet at 13,000 rpm for 10 min at 4°C. Transfer supernatant (This is the Cell Lysate) to a new tube.
5. Assay protein concentration of the Cell Lysate. The Cell Lysate can be used immediately or freeze at -80°C for future use.

B. “Pull Down” JNK Using c-Jun Fusion Protein:

6. For each assay, add 20 μ l c-Jun Fusion Protein Beads to 200 μ l Cell Lysate (~ 50 - 400 μ g total protein). Incubate with gentle rocking overnight at 4°C.
7. Microcentrifuge (14,000 rpm) for 30 sec at 4°C. Remove Supernatant. Wash pellet twice with 0.5 ml of 1X Kinase Extraction Buffer and one time with 0.5 ml Kinase Assay Buffer. Keep on ice.

C. Kinase Assay:

8. Suspend pellet in 50 μ l Kinase Assay Buffer. Add 1 μ l of 10 mM ATP. Incubate for 30 min at 30°C.
9. Add 30 μ l 3X SDS-PAGE Buffer (not provided).
10. Boil the samples for 3 min. Microcentrifuge for 2 min.
11. Load the supernatant (20 μ l) on 12% SDS-PAGE. Alternatively, the supernatant may be stored at -20°C for future use.

D. Western Immunoblotting:



12. Perform Western blot analysis using the rabbit anti-Phospho-cJun (Ser 73) Specific Antibody at 1:1000 dilution. A 35 kDa band corresponding to the phosphorylated c-Jun protein should be detected in JNK activated samples.

For research use only! Not to be used in humans.

Our promise

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For more details, please visit <http://www.apexbt.com/> or contact our technical team.

Tel: +1-(832)696-8203

Fax: +1-832-641-3177

Email: sales@apexbt.com